I Istruzioni per l'uso - **Logic Sphero**
GB Instructions for use - **Logic Sphero**
D Gebrauchsanleitung - **Logic Sphero**
E Instrucciones de uso - **Logic Sphero**
CZ Návod k použití - **Logic Sphero**

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English
GB
GENERAL WARNING

- Dental Tech implants are only to be used by specialised medical personnel who are properly trained.
- The use of components which are not part of the system produced by Dental Tech may compromise the functionality of the implants.
- The implants need to be stored in their original packaging in a dry place out of direct sunlight; before inserting them check the packaging and the implant to make sure they are undamaged.
- Do not use implants which are not sterile or which have been previously used. Dental Tech shall not be held liable for implants which are re-sterilised by the user and shall not accept returns of implants contained in original packaging which has been tampered with.
- Do not use the implant after the expiry date shown on the package.
- The information on the use of Dental Tech implants is indicated in the surgical protocol and product catalogue and reflect the current state of progress known to Dental Tech at the time the product is put on the market. This information must only be considered as a recommendation and does not exonerate the user from ensuring that the product is adequate for the intended use.
- Processing and application of the product, occurs outside of our control and liability. We shall not be held liable for any damage resulting from this.
- Dental Tech shall not be held liable for failure to comply with the above rules.

INDICATIONS

- The Logic Sphero implants are made of Titanium alloy (Ti-6Al-4V) and must only be used when the dimensions of the ridge are not sufficient to allow the use of conventional implants with a greater diameter. Moreover, they must only be used to stabilise total mobile prostheses in the inferior arch. Therefore at least 4 Logic Sphero implants needs to be inserted in the intraforaminal region which are parallel and equidistant from each other, to obtain a stable connection of the prosthesis. Any inclination of each single implant, compared to the parallel axis, must not exceed 8°. There must also be a minimum distance of not less than 7 mm between the implant sites.

Important: Instructions to follow for inserting implants.

- To prepare the implant site use Ø2.0 osteotribes or shafts of the implant diameter.
- Constantly spray the osteotribes and shafts with plenty of sterile physiological solution which has been previously cooled (temperature between 4° and 10°).
- It is important to note that the two depth marks, corresponding to the two available endosseous heights, are indicative only if a gum flap has been made for a total preparation of the site height (the selected height notch must reach the cortex). In these cases it is necessary to wait for healing of the tissue before applying the prosthetic load.
- For manual insertion the device must be connected to the contra-angle by connection for contra-angle (RDS225). In this case do not exceed 25 rpm and 55 Ncm torque.
- Couple the Dental Tech prosthetic system composed on the titanium crown POR225 with the silicone o-ring ORG225, or alternatively to specific retentive caps and related steel cases.
- It is advisable to replace the o-ring after evaluation of the prosthesis hold. If retentive caps are used the manufacturer's instructions should be followed.
- It is possible to plan the clinical case by evaluating the endosseous height of the implants to insert by using x-rays and specific display device; however, this instrument is not to be used for taking measurements, but must only provide indications on the form of the implants.
- To take impressions use transfers and relative analogues.

CONTRAINDICATIONS

- General serious health problems of the patient;
- Diabetes;
- Psychosis, immunodepression, empathies;
- Heart or vascular disease;
- Patient with poor oral hygiene habits;
- Insufficient or poor quality bony tissue (which could compromise the stability of the implant);
- Drug, alcohol or tobacco abuse;

- Chronic mandibular osteitis;
- Chronic or acute infectious diseases;
- Localised bone tumours or irradiation of the implant site;
- Systemic diseases;
- Bruxism;
- Pregnancy;
- Advanced osteoporosis;
- Current treatment with anticoagulants or immunosuppressants.

SIDE EFFECTS

To minimise any side effects the patient must refrain from heavy physical activity after the operation for at least 4 days. Temporary side effects may occur, including: local pain, swelling and local inflammation. Long term side effects are also possible, including: 1) chronic pain connected to the implant; 2) temporary or permanent damage to nerve vascular bundles with paresthesia or anaesthesia; 3) bone absorption at the mandibular maxillary torus; 4) maxillary bone fractures, nerve lesions.

PACKAGING AND STERILITY

The Dental Tech implant is supplied in a sterile package which, if undamaged, guarantees the protection of the implant against external agents, and if correctly stored, its sterility.

The package supplied by Dental Tech is composed of the following items:

- Cardboard box with enclosed instruction sheet, removable label to apply to the patient's chart and information sheet for the patient;
- Blister in plastic labelled and sealed to guarantee the sterility of the implant;
- Plastic ampoule (which contains the implant and mounting device).

REMOVING THE IMPLANT FROM THE STERILE PACKAGE

To remove the implant from the sterile package it is necessary to follow aseptic and sterility rules and work in a sterile environment. The package must only be opened at the time of the operation and after completing the preparation of the implantation site. Avoid contact with non-sterile items during all the phases. To remove the implant proceed as follows: 1) open the blister; 2) let the ampoule fall onto a sterile cloth; 3) open the ampoule by removing the cap with connected implant; 4) insert the implant in the implantation site using the ampoule cap as described in the surgical protocol.

DOCUMENTATION AND IDENTIFICATION

Dental Tech recommends complete clinical, x-ray and statistical documentation. The item can be identified based on the item code and batch number shown on the package.

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