Via G. Di Vittorio, 10/12
20826 Misinto (MB) - Italy

Tel. + 39 02 967.20.218
fax +39 02 967.21.269

info@dental-tech.it
www.dentaltechitalia.com

Italiano

I

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GB

INTENDED USE / INTENDED FUNCTION

Dental implants are intended to replace missing or corrupted teeth, • that are not possible to be repaired, replaced or compensated by other means; • where other solutions have an undesired impact to sound teeth, or • where implants are desired for obtaining an optimal cosmetic result. In general, COMPRESSIVE implants are intended for surgical placement in the upper or lower jaw to provide an anchorage for prosthetic superstructures for tooth restorations or as a terminal or intermediary abutment for fixed or removable bridgework, and to retain overdentures.

INDICATIONS

The medical indications for the use of a COMPRESSIVE implant are:

- Loss of teeth / missing teeth,
- replacement of damaged or ill teeth.

The concrete disease, injury, physiological condition or traumatic event leading to the loss of a tooth or to the necessity of tooth removal are manifold and do not matter, as long they are not explicitly listed in the contraindications.

RANGE OF APPLICATION

The COMPRESSIVE implant is a one-piece implant with compressive threads. It is used for multiple unit restorations with immediate loading in the upper and lower jaws with adequate bone tissue. Abutment direction can be adjusted up to 15° relative to the implant axis. Can be used with caution to create single restorations in situations where high primary stability is achieved on placement.

CONTRAINDICATIONS

The contraindication to place implants will be evaluated by professionals on a case-by-case basis, with the greatest caution. Preoperative diagnosis is necessary to identify threats to the patient, related to the procedure of the implant placement, as well as factors that may affect the possibilities of healing of the bone and surrounding soft tissues. Contraindications can be separated into absolute and relative contraindications:

ABSOLUTE CONTRAINDICATIONS :

- Myocardial infarction: within six months of an attack,
- Cerebral infarction and cerebral apoplexy: In cases where the condition of the disease is serious and the patient is concurrently taking anticoagulants,
- Severe immunodeficiency,
- Patients who are undergoing strong chemotherapy,
- Severe neuropsychiatric disease, mental disability, and narcotic drug addicts,
- Patients who are concurrently taking bisphosphonates,
- Allergies or hypersensitivities to chemical ingredients of material used: titanium alloy.

RELATIVE CONTRAINDICATIONS:

- Diabetes (particularly insulin-dependent),
- Angina pectoris (angina),
- Seropositivity (absolute contraindication for clinical AIDS),
- Significant consumption of tobacco,
- Certain mental diseases,
- Radiotherapy to the neck or face (depending on the zone, quantity of radiation, localization of the cancerous lesion etc.),
- Certain auto-immunes diseases,
- Drug and alcohol dependency,
- Pregnancy,
- Certain diseases of the mucous membranes of the mouth,
- Bruxism,
- Periodontal diseases (loosening of the teeth); it is necessary to clean up the gums and stabilize the disease first,
- An unbalanced relationship between the upper and lower teeth,
- Poor hygiene of the mouth and teeth,
- An insufficient quantity of bone,
- Infections in the neighboring teeth (pockets, cysts, granulomas), major sinusitis,

SIDE EFFECTS, COMPLICATIONS WITH COMPRESSIVE IMPLANTS

Immediately after the insertion of dental implants, activities that demand considerable physical exertion should be avoided. Possible complications following the

insertion of dental implants are:

Temporary symptoms: pain, swelling, phonetic difficulty and gingival inflammation.

More persistent symptoms: chronic pain in connection with implants, permanent paraesthesia, dysesthesia, loss of maxillary / mandibular ridge bone, localized or systemic infection, oroantral or oronasal fistula, unfavorably affected adjacent teeth, fracture of implant, jaw, bone or prosthesis, aesthetic problems, nerve damage, exfoliation, hyperplasia.

WARNING

Products must be secured against aspiration when handled intraorally. Aspiration of products may lead to infection or unplanned physical injury. If you want to protect it, use rubber dam!

Do not exceed recommended insertion torque, as it might cause bone necrosis or system components fracture.

Beside the mandatory precautions for any surgery such as of asepsis, during drilling in the jaw bone, one must avoid damage the nerves and vessels by referring to anatomical knowledge and preoperative medical imaging (e.g. radiographs).

Failure to recognize actual lengths of drills relative to radiographic measurements can result in permanent injury to nerves and other vital structures. Drilling beyond the depth intended for lower jaw surgery may potentially result in permanent numbness to the lower lip and chin or lead to hemorrhage in the floor of the mouths.

Do not use device if the primary package has been damaged or previously opened.

Do not use damaged or blunt instruments for implantation.

CAUTIONS / PRECAUTIONS

- Sterile handling is essential. Never use potentially contaminated components. Contamination may lead to infection.
- Do not resterilize Dental Implants.
- Avoid any contact of the implant with foreign substances prior to their use. Do not touch the endoseal part of the implant.
- Implants should be used according to their expiration date.
- If implants not assembled any more with holder and just moving into the blister, DO NOT USE this implant because surface already contaminated by plastic particles.
- Treatment by means of implants may lead to loss of bone, biologic and mechanical failures, including fatigue fracture of implants.

CLEANING AND DISINFECTION

Implants are delivered sterile and for single use only prior to the labeled expiration date. They must not be cleaned and sterilized.

STERILIZATION

Implants are delivered sterile. The intact sterile packaging protects the sterilized implant from external influences and if stored correctly, the packaging ensures sterility up to expiration date.

The sterile packaging must not be opened until immediately prior to insertion of the implant. When removing the implant from sterile packaging, rules of asepsis must be observed. The sterile packaging must not be opened until immediately prior to insertion of the implant.

INSERTION OF THE IMPLANT

The implant is removed from the sterile packaging immediately prior to the introduction and shall be stably installed in the implant bed prepared. Be sure to install it securely immediately. COMPRESSIVE Implant can be placed either manually with the ratchet or with the aid of the handpiece.

DOCUMENTATION AND IDENTIFICATION

Dental Tech recommends a complete clinical radiographic and statistical documentation. The item can be identified by its item code, its batch number, both on the packaging

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