- I** Foglio istruzioni - **Pilastro MUA**
- GB** Information leaflet - **MUA Abutment**
- D** Gebrauchsanweisung - **Abutment MUA**
- E** Hoja de instrucciones - **Pilar MUA**
- CZ** Návod k použití - **Abutment MUA**

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GB

### GENERAL WARNING

- Dental Tech abutments MUA and healing cuffs shall be used by specialised medical staff and dental technicians only.
- The use of prosthesis components not belonging to the endosseous system manufactured by Dental Tech, may impair its functionality.
- The product is manufactured and applied out of our control and responsibility. The manufacturing company cannot be held responsible for any deriving damage.
- Prior to using any Dental Tech product inspect the packaging and device for damages.
- At the end of its useful life the medical device must be disposed of in accordance with national laws regarding waste disposal.
- Reuse of the device is prohibited. If the operator decides to violate this prohibition, the patient may be exposed to risks of infections, diseases transmitted through contact with blood or bodily fluids, and deformation of the device, make it impossible to match with related medical devices (implants, instruments, etc.)
- Dental Tech cannot be held responsible for any damage resulting from inobservance of this warning.

### INDICATION

#### - FOR MUA ABUTMENT

- MUA Abutment for multiple elements is suitable to manufacture the following:
  - Screwed bridges
  - Toronto/Branemark® screwed bars
  - Ackermann screwed bars
  - Dolder screwed bars
  - Drilled screwed bars (Mesostructured) also suitable for removable prostheses supported by bolts
- Supports of temporary and permanent structures for immediate load

#### - FOR HEALING ABUTMENTS

- Dental Tech healing abutments shall be used exclusively matching MUA Abutment
- Healing cuffs are used to make an ideal profile to the secondary part and are suitable to heal soft tissues at the second surgical stage.
- Leave them in the oral cavity for 10/15 days.
- To make screwing/unscrewing easier, use the short screw-driver code GMM250.

#### - IMMEDIATE LOAD

In the event of immediate load on Abutment MUA, connect it to implants just after its insertion, stitch soft tissues around the abutment where you can make the impression by using the IMPLANT IMPRESSION specifically designed for abutment MUA suitable for the corresponding prosthetic platform.

#### - SPREAD LOAD

In the event of spread load, connect the abutment MUA to implants at the second surgical stage, stitch soft tissues around the abutment where the special HEALING CUFF, suitable for the corresponding prosthetic

platform, is then screwed until soft tissues heal. Once healed, make the impression by using the IMPLANT IMPRESSION specifically designed for abutment MUA suitable for the corresponding prosthetic platform.

#### - DEVELOPING THE IMPRESSION

In case of immediate load or in case of spread load, the technician shall develop the impression and make the working sample by using PLASTER ANALOGS specifically designed for abutment MUA suitable for the corresponding prosthetic platform.

#### - MAKING PROSTHETIC STRUCTURE

The technician shall make the prosthetic structure by using the PLASTIC SLEEVE or TITANIUM ABUTMENT specifically designed for abutment MUA suitable for the corresponding prosthetic platform, or also with CAD-CAM prosthetic techniques, making use of digital libraries.

### SPECIFIC WARNING

- Screw the abutment MUA to the implant by using the torque wrench, tightened at 30 Ncm as well as the special adapter code TW0080.
- Do not screw the abutment MUA to the implant by using unsuitable tools.

### CONTRAINDICATIONS

#### - FOR MUA ABUTMENT

The abutment MUA is not suitable for cementum-covered prostheses, as well as for manufacturing screwed single elements since it does not have any anti-rotational device.

#### - FOR HEALING CUFFS

- Do not screw or unscrew with unsuitable tools.
- Use it for the above-mentioned purpose only.
- Do not use for soldering, welding or permanent assemblages.
- Do not modify or drill.

### PACKAGING AND STERILITY

- This device is disposable.
- "Autoclave moist heat sterilisation must be carried out by using an validated and routinely checked method pursuant the requirements as per ISO 17665-1:2007 "Sterilization of health care products. Requirements for validation and routine control of moist heat sterilization in health care facilities". An example of cycle may be represented by autoclave sterilisation able to get rid of any air pocket within the sterilisation chamber for 20 minutes at a temperature equal to 134°C and pressure equal to 2 atm.".

### WARNING

**Pouches unsuitable for autoclave sterilisation.**

*N.B. Cleaning and sterilisation shall be carried out by qualified personnel, under the supervision of a doctor and by using properly calibrated and faultless machines.*

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