Italiano

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English

GENERAL WARNING

- Dental Tech abutments shall be used by specialised medical staff and dental technicians only.
- The use of prosthesis components not belonging to the endosseous system manufactured by Dental Tech, may impair its functionality.
- The product is manufactured and applied out of our control and responsibility. The manufacturing company cannot be held responsible for any deriving damage.
- Prior to using any Dental Tech product inspect the packaging and device for damages.
- At the end of its useful life the medical device must be disposed of in accordance with national laws regarding waste disposal.
- Reuse of the device is prohibited. If the operator decides to violate this prohibition, the patient may be exposed to risks of infections, diseases transmitted through contact with blood or bodily fluids, and deformation of the device, make it impossible to match with related medical devices (implants, instruments, etc.)

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GB Information leaflet - Temporary abutment

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- Dental Tech cannot be held responsible for any damage resulting from inobservance of this warning.

DIRECTIONS

- Temporary abutments MKS - HKS are suitable for permanent cementum-covered prostheses, prep temporary abutments PKE - PKR - PPE - PPR- HPP - HPC - HKE - HKR for cementum-covered and taper prostheses.
- Temporary abutments can be adjusted within the oral cavity and, thanks to its anti-rotational system you can use the indirect technique by using impressions keeping the surgical position of the implant in the impression, and analogs reflecting the inside part and the neck of the implant.
- In the event of any adjustment in the oral cavity, sprinkle generously. To adjust temporary abutments we recommend using drills specifically designed for titanium.
- Temporary abutments can be adjusted to reach up to a minimum thickness of 0.5 mm in the coronal portion and 0.6 mm in the cervical portion.
- After manufacturing the prosthesis, before permanent and/or temporary casting, we suggest protecting the abutment screw with special material.
- Soft tissues needs pre-conditioning through special healing abutments.
- Tighten abutment screws with Dental Tech torque wrench at 20 N cm and the special adapter for torque wrench code TW0001 C - TW0001 L.
- Do not weld or solder.
- If the transmucosal part is to be adjusted, do not exceed the connecting area to the implant.

CONTRAINDICATIONS

- Do not screw or unscrew the abutment screw with unsuitable tools.
- Use for the above-mentioned purpose only.

PACKAGING AND STERILITY

- This device is disposable.
- "Autoclave moist heat sterilisation must be carried out by using an validated and routinely checked method pursuant the requirements as per ISO 17665-1:2007 "Sterilization of health care products". Requirements for validation and routine control of moist heat sterilization in health care facilities". An example of cycle may be represented by autoclave sterilisation able to get rid of any air pocket within the sterilisation chamber for 20 minutes at a temperature equal to 134°C and pressure equal to 2 atm".

WARNING

Pouches unsuitable for autoclave sterilisation

N.B. Cleaning and sterilisation shall be carried out by qualified personnel, under the supervision of a doctor and by using properly calibrated and faultless machines.

THE PEEK PROVISIONAL ABUTMENTS HAVEN'T TO REMAIN INSIDE THE MOUTH MORE THEN 180 DAYS.

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