## Italiano

## I

### USO PREVISTO

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## GB

### INTENDED USE

Cobalt chrome overcast abutments are used as secondary components connected to dental implants, they are used as prosthetic components in implant-prosthetic rehabilitations of clinical cases previously treated with the placement of endo-bone implants. They act as a mechanical connection between the endosseous dental implant and the dental crown.

### INDICATION FOR USE

Cobalt chrome overcast abutments are secondary components used to produce prosthetic artefacts (cemented and/or screwing) on implants, as an alternative to the other prosthetic components available, when the characteristics of dental restoration must be particularly adapted to the aesthetic and functional specifications of the treated clinical case. Cobalt chrome abutments are designed by the customer through traditional processing, or through an intraoral scan and designing the shape using dental software. Cobalt chrome abutments stumps can be coated directly in ceramic. They are also used to make screwed bars, Toronto bars or Toronto Branemark® and their purpose is to ensure precise coupling with the dental implant.

### WARNINGS

- Dental Tech cobalt chrome overcast abutment should only be used by specialized dentist and specialized dental technicians.
- The devices must be used for the sole purpose of completing an implant-prosthetic treatment carried out with endosseous implants with a corresponding connection (see compatibility table).
- The use of prosthetic components not belonging to the endosseous system produced by Dental Tech, could impair its functionality.

- The processing and application of the products is outside our control and responsibility. Liability in the event of damage resulting from them is excluded.
- Check the integrity of the device before use. Do not use the product if the packaging is damaged or opened.
- At the end of its shelf-life, the medical device must be disposed of in accordance with the procedures laid down by the national waste laws in force.
- The reuse of the device is prohibited: if the operator decides to break this prohibition, it may expose the patient to risks of infection, diseases caused by contact with the patient's blood or bodily fluids, deformation of the device and consequent impossibility of coupling with associated medical devices (dental implants, abutments, etc.).
- Overcast with alloys suitable for the overcast of metals for dental use in cobalt chrome that have a melting interval of 1075-1150 °C, without exceeding the deformation temperature between 1120-1200 °C, in order not to alter the characteristics of the preformed base.
- Tighten the abutments exclusively with the dedicated fixing screw and by dynamometer tool at 20Ncm.
- Dental Tech disclaims all responsibility in case of non-compliance with the above-mentioned warnings.

### CONTRA-INDICATIONS

- Do not use in case of known nickel allergy.
- Do not screw or unscrew the fastening screw or screwable pillars with occasional tools.
- Do not use melting alloys that may compromise their biocompatibility.
- Cobalt chrome overcast abutment, with Dental Tech fastening screw is not suitable for the realization of immediate temporary devices.
- Do not use the overcast abutment in the patient's oral cavity before it has been transformed and modified.

### SIDE EFFECTS

Temporary side effects of prosthetic implantology interventions can be pain, swelling, hematoma, pronunciation problems, and soft tissue inflammation. The success of implant-prosthetic therapy is closely linked to the correctness of all the stages leading to its completion, including the planning and implementation of the prosthesis. The correct and careful use of prosthetic components is a very important aspect that favors the reduction of the risks of failure of implant-prosthetic treatment, as well as reducing the extent of temporary side effects. It is required to report to Dental Tech S.r.l. any cases of failure.

### PACKAGING

- The cobalt chrome overcast abutment is supplied, non-sterile, in single bags with its fixing screw and overcast component.
- The device is disposable.
- The device after being transformed and before being inserted into the oral cavity, it must be cleaned, disinfected and sterilized (along with the fixing screw).
- Packaging bags not suitable for sterilization in autoclave.

### CLEANING, DISINFECTION AND STERILIZATION

- Cleaning, disinfection and sterilization operations must be carried out by qualified personnel under the direct responsibility of the doctor; sterilization must be carried out through the use of regularly calibrated machinery and in perfect state.
- Wear disposable gloves during all cleaning, disinfection and sterilization operations.

Dental Tech has validated the manual cleaning and disinfection method below. Do not use automated cleaning and disinfection methods because they are not validated.

### MANUAL CLEANING

- Before use of Dental Tech abutment, place the abutment into ultrasound bath with a peracetic acid based solution at concentration of 2% (NO GLUTARALDEHYDE OR SODIUM HYPOCHLORITE), as long as 18 minutes.
- Afterward rinse carefully with drinking water.

### MANUAL DISINFECTION

- Place the abutment into a container with a peracetic acid based solution at concentration of 4% (NO GLUTARALDEHYDE OR SODIUM HYPOCHLORITE), as long as 15 minutes.
- Rinse generously with drinking water.
- Examine the abutment and make sure there are no organic remains. Carefully scrub the outer parts with a non-metal rigid bristled brush.

### MANUAL RINSE

- Place the abutment into ultrasound bath with distilled water, and wash it for approx. 18 minute and then rinse carefully with distilled water.

### DRY

- Dry perfectly in aseptic conditions (gauze, paper, filtered air...) and store the abutments individually bagged using Medical Paper C and FILM PET/PP Peel (ISO Reference 11607-1).

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### STERILIZATION (according to ISO 17665-1)

Dental Tech validated the following Autoclave moist heat sterilization cycle:

- Pre-vacuum time according to autoclave builder specifications.
- Temperature: 134 °C.
- Holding time: 3 minutes.
- Drying time: minimum 10 minutes.

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