Istruzioni per l'uso - Impianti dentali
Instructions for use - Dental implants
Gebrauchsanleitung - Implantate
Instrucciones para uso de - Implantes dentales
Návod k použití - Zubní implantáty

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English GB

IMPLASSIC - IMPLOGIC® implants are made of grade 5 titanium with BWS sandblasted surface (aluminium oxide treatment and subsequent treatment with acids).

CAUTIONS

- Dental Tech endosseous implants shall be used by specialised and properly trained medical staff only.
- The use of prosthesis components not belonging to the endosseous system manufactured by Dental Tech, may impair its functionality.
- The information on how to use Dental Tech implants are described in the implant protocol as well as in the catalogue of products and they mirror the state of the art procedure known to Dental Tech when the product has been marketed. This information is meant only as recommendations and does not prevent the user from making sure that the product is suitable for its intended purpose.
- The product is manufactured and applied out of our control and responsibility. The manufacturing company cannot be held responsible for any deriving damage.
- Do not insert non-sterile or previously used implants. The manufacturing company cannot be held responsible for Dental Tech implants which have been sterilised again by the user.
- Do not use the implant after the expiration date on the packaging.
- At the end of its useful life the medical device must be disposed of in accordance with national laws regarding waste disposal.
- Reuse of the device is prohibited. If the operator decides to violate this prohibition, the patient may be exposed to risks of infections, diseases transmitted through contact with blood or bodily fluids, and deformation of the device, make it impossible to match with related medical devices (implants, instruments, etc.)
- Dental Tech cannot be held responsible for any damage resulting from inobservance of this warning.
- Implants shall be kept in their original packaging in a dry place and out of the sunlight; prior to its insertion make sure either the packaging or the implant is not damaged.

DIRECTIONS

Dental Tech endosseous system are intended to be used when treating patients suffering from partial or full edentulism on both arches (upper and lower). Dental Tech implants can be turned into prostheses through single crowns, bridges, partial prostheses, immediate or total overdentures by using secondary components Laser marking and the colour code on secondary components or by the laboratory enables them to be quickly matched to the implant

CONTRAINDICATIONS

- Serious health problems in the patient;
- Diabetes;
- Psychosis, immune-depression;
- Cardiovascular disorders
- Poor oral health in the patient
- Insufficient or poor-quality hard tissue (which may impair the stability of the implant)
- Drug, alcohol and tobacco abuse ;
- Severe maxillary osteitis
- Severe or acute contagious diseases;
- Localised bone cancer or irradiation of the implant site;
- Systematic diseases;
- Pregnancy ;
- Advanced osteoporosis;
- Ongoing treatments with anticoagulants or immunosuppressive drugs;
- Haemopathy.

SIDE EFFECTS

To minimise side effects, the patient shall not do any strenuous physical activity just after the operation. The patient may experience temporary side effects such as: local pain, swelling, local inflammation. The patient may

also experience long-term side effects such as: 1) severe pain connected to the implant; 2) temporary or permanent damage to nervous bundles with paraesthesia or anaesthesia; 3) bone recession at the jawbone crest; 4) jawbone fractures, nerve injuries.

PACKAGING AND STERILITY

Dental Tech endosseous implants are supplied in sterile packaging which enables - if undamaged - to protect the implant from external agents and - if properly stored - its sterility. The packaging provided by Dental Tech is made up of the following parts:
• Cardboard box containing a leaflet, removable label to apply onto the patient's file, the patient information leaflet;
• Plastic blister labelled and sealed assuring the sterility of the implant;
• Plastic ampoule (containing the implant, assembling device and the cover screw)

HOW TO REMOVE THE IMPLANT FROM THE STERILE PACKAGING


To remove the implant from the sterile packaging you need to meet the sterility and asepsis requirements and work in a sterile environment. The packaging shall be opened only when surgically operating after fully preparing the receiving site. At all stages avoid any contact with non-sterile elements.

To remove the **ImpLassic HX3** implant from its packaging for insertion, proceed as follows:

- 1) open the blister;
- 2) remove the ampoule, keeping the cap (A) pointed upwards Fig.1;
- 3) remove the cap (A) and reveal the position for the implant;
- 4) attach the contra-angle screwdriver (HHP2521 - HHP2525 - HHP2529) on the contra-angle, to remove the implant from the ampoule;
- 5) insert and screw the implant in the alveolus until it is completely inserted.
Warning Max. 50 Ncm torque during insertion.

To remove the **ImpLassic FT3 / ImpLassic FTP / ImpLogic® AT / ImpLogic® GII** implant from its packaging for insertion, proceed as follows:

- 1) open the blister;
- 2) remove the ampoule, keeping the cap (A) pointed upwards Fig.1;
- 3) remove the cap (A) and reveal the position for the implant;
- 4) attach the contra-angle screwdriver (ICO325-ICO3727-ICO3732 and ICOFTP100-ICOFTP200-ICOFTP300 for ImpLassic FTP) on the contra-angle, to remove the implant from the ampoule;
- 5) insert and screw the implant in the alveolus until it is completely inserted.
Warning Max. 70 Ncm torque during insertion.

To remove the closure screw from the packaging, proceed as follows: remove the closure screw from the ampoule, keeping it pointing up towards (B),

using the hexagonal screwdrivers or the contra-angle with connector (GCG0024 - GCG0030). Fig.1
To remove the **ImpLassic TR2** and **TRW2** implant from its packaging for insertion, proceed as follows:

- 1) open the blister;
- 2) remove the ampoule, keeping the cap (A) pointed upwards Fig.1;
- 3) remove the cap (A) and reveal the position for the implant;
- 4) attach the contra-angle screwdriver (CML037-CML037) to the contra-angle, to remove the implant connected to the mount, from the ampoule.
- 5) insert and screw the implant in the alveolus until it is completely inserted.
Warning Max. 70 Ncm torque during insertion.

To remove the closure screw from the packaging, proceed as follows: remove the closure screw from the ampoule, keeping it pointing up towards (B), using the hexagonal screwdrivers or the contra-angle with connector (GCG0024 - GCG0030). Fig.1

Consult the dedicated SURGICAL MANUALS for further clarifications.

DOCUMENTATION AND IDENTIFICATION

Dental Tech recommends a complete clinical radiographic and statistical documentation. The item can be identified by its item code, its batch number, both on the packaging.

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