



Dental Tech: a long road dedicated to implantology

Dental Tech is a long-standing and responsible company, in constant growth both in Italy and abroad, located in Misinto (MB). With 40 years of experience manufacturing implantable medical devices, Dental Tech has become a brand synonymous with quality and reliability. Dental Tech aims to provide a safe and versatile line of implants systems and to ensure that the wide range of products meets the needs of professionals and to cope with the different clinical cases, always maintaining a high standard of quality and the right price. With the aim to set the standards of tomorrow and improving existing ones, Dental Tech has a technical staff constantly employed in research and development to offer the best solutions to professionals and patients. The result is a constant know-how, always in step with the times to be the point of reference, able to satisfy the needs of the dentist and the well-being of the patient.

Dental Tech's strength is a team of experts and highly professional sales force that support clinicians with their daily requirements. Besides, thanks to the synergy with qualified clinical dentists are trained in the use of the implant system, through advanced courses suitable for the implementation of the clinical performance of the participants.

The wide Dental Tech's implant line offers conical and cylindrical systems, internal and external connection, transmucosal and submerged systems, short implants and one-component mini implants (Logic Sphero). The implants are available (depending on the implant line) in lengths from 6 to 16 mm and in different diameters (Ø3.25, Ø3.75, Ø4.25, Ø4.75 and Ø5.5).

The implant line is subjected to a double process of sandblasting and acid etching that makes it possible to obtain optimal values of roughness creating the strongest fibrin adhesion to the



surface and facilitating the bone healing process by significantly reducing the time.

After the surface treatment and the classic washings, Dental Tech implants are additionally cleaned with Argon Cold Plasma to minimize carbon contamination subsequently and minute controls are performed on the fixture with scanning electron microscopes (SEM). Packaging in controlled environments (clean room ISO 7), packaging impermeable to micro-organisms and gamma ray sterilisation process guarantee the creation of products that are extremely safe for users and their patients.

Dental Tech operates in compliance with the rigorous requirements of European Directive 93/42/EEC, is certified by TUV Product service CE 0123 for medical products, and our quality management system is ISO certified in accordance with UNI EN ISO 9001:2008 EN ISO 13485:2012+AC:2012.



Contact

Address: Via G. Di Vittorio, 10/12
Misinto (MB) 20826

Phone: +390296720218 / +390296720174 / +390296720178

Fax: +39 02 96721269

info@dental-tech.it

www.dentaltechitalia.com