

• Dental Tech

Dental Tech aims to provide a safe and versatile line of implants systems and to ensure that the wide range of products meets the needs of professionals.

The new ImpLogic® AT (Active Thread) is a tapered implant that, thanks to its special spiral design, facilitates the users in the realization of Ridge expansion procedures. The exceptional self-tapping power of the thread provides an excellent bone condensation and a high primary stability even in very complex clinical cases. ImpLogic® AT is recommended in cases of post extraction implants and in case of poor quality bone.

The main characteristics of the implant system are:

- **SPIRAL DESIGN:** the unusual spiral design simplifies the procedures of Ridge expansion.
- **RISK REDUCTION:** less risk of damaging adjacent teeth and perforation of the lingual and/or buccal cortical plates.
- **BONE MAINTENANCE OVER TIME:** allows a greater reduction of bone osteotomy, which results in lower bone loss and reduced surgical trauma.
- **SELF-TAPPING COIL:** exceptional self-tapping capability, which provides improved bone condensation and increased primary stability, even in highly complex clinical cases.
- **OPTIMAL CHOICE OF POSITIONING:** allows a change in direction in order to achieve the optimum position of restoration, especially in post-extraction sites.

ImpLogic® AT is available in diameters \varnothing 3.75, \varnothing 4.25 and \varnothing 4.75 and in lengths 8, 10, 11.5, 13 and 16 mm*.

*No for \varnothing 4.75

www.dentaltechitalia.com
info@dental-tech.it

