

ImplantBook®

2019 | Global guide for Dealers and Dentists



DENTAL IMPLANTS

AROUND IMPLANTOLOGY

CONGRESSES



• **Dental Tech**

ImpLassic HX3



Manufacturer: ImpLassic HX3, produced by **Dental Tech**, replaces and updates the previous implant line with external connection ImpLassic HX2.

Morphology: ImpLassic HX3 is a cylindrical body implant and conical apex. Modulating the surgical procedure is indicated in all types of bone and even in the case of non-compact bone is able to achieve excellent primary stability. The external hexagon connection makes it particularly suitable for the realization of Full Arch prostheses on Fixtures and screwed prostheses in general.

Surface: The BWS® surface, made by Dental Tech through a sandblasting and etching process, provides an extremely clean surface and an optimum roughness that promotes cell adhesion.

Available Implant Lengths (mm): 8* / 10 / 11.5 / 13 / 16**
* No for Ø3.25 ** Only for Ø4.0

Available Implant Diameters (mm): 3.25 / 4.0 / 4.5

Seating Position To The Bone Crest: Crestal level.

Abutment-Implant Connection Type:

Universal connection with external hexagon.
Platform 3.5 with a hexagon of 2.4 for Ø3.25.
Platform 4.1 with a hexagon of 2.7 for Ø4.0 and Ø5.0.

Index type abutment repositioning: The implant system's index for the abutment repositioning is hexagonal.

info@dental-tech.it // www.dentaltechitalia.com

• **Dental Tech**

ImpLassic FTP



Manufacturer: ImpLassic FTP, produced by Dental Tech, completes the safe and wide range of products that the company makes available to meet the needs of all professionals.

Morphology: ImpLassic FTP is a conical implant system with spiral profile with hybrid progress: flat and radiating towards the root, triangular-shaped externally, for greater penetration into incompletely prepared sites. Modulating the surgical procedure is indicated in all types of bone and even in the case of non-compact bone is able to achieve excellent primary stability.

Surface: The BWS® surface, made by Dental Tech through a sandblasting and etching process, provides an extremely clean surface and an optimum roughness that promotes cell adhesion.

Available Implant Lengths (mm): 8* / 10 / 11.5 / 13 / 16**
* No for Ø3.75 ** Only for Ø4.75

Available Implant Diameters (mm): 3.75 / 4.25 / 4.75

Seating Position To The Bone Crest: It's a submerged implant with a 0.5 / 1 mm subcrestal placement

Abutment-Implant Connection Type: Implant system with 6° conometric connection to facilitate positioning of the abutment.

Index type abutment repositioning: The implant system's index for the abutment repositioning is hexagonal.

info@dental-tech.it // www.dentaltechitalia.com




DentalTech[®]

DENTAL IMPLANT SYSTEM

DESIGNED & MADE

in Italy



DENTAL TECH SRL
Via G. Di Vittorio, 10/12
20826 Misinto (MB), Italy

Phone +39 02 967.20.218
info@dental-tech.it
www.dentaltechitalia.com

