

ImplantBook®

2021 | Global guide for Dealers and Dentists



Compendium

DENTAL IMPLANTS

AROUND IMPLANTOLOGY

IMPLANTBOOK REPORT



www.Implant-Book.com
ImplantBook@Infodent.com



DentalTech®

DENTAL IMPLANT SYSTEM

MADE IN ITALY



DENTAL TECH IMPLANT LINE

Safe, versatile and for all clinical needs

+39 02 96720174

info@dental-tech.it

www.dentaltechitalia.com

Transmucosal Abutment

Shape Tissue Effect



Dental Tech has expanded its wide range of products with a new prosthetic component that promotes the conditioning of peri-implant tissues. It is applied at the same time as the surgical phase of implant placement, maintaining freedom of flexibility during the prosthetic phases. The transmucosal abutment radically simplifies the prosthetic procedure, since it moves the engagement platform of the external connection implants from the bone level to the soft tissue level. Moreover it simplifies the impression taking since the connection is clearly visible.

It is prothesisable both in immediate loading, in a variety of secondary components suitable for the planned prosthetic project, and in deferred loading, preserved by a protective cap that leaves the soft tissue undisturbed to optimize healing and complete Implant integration. The abutment's external connection platform offers a wide range of prosthetic solutions to better manage the prosthetic needs of each project. Once positioned according to the prosthetic emergency, it is no longer removed, guaranteeing the stability of the surrounding peri-implant tissues. Acting on the transmucosal abutment with a Torque of 30Ncm, perfect matching between the abutment-implant interface is ensured, thus preventing the migration and bacterial colonization on the "pumping effect" action.

The abutment is available for FTP, ImpLogic® and ImpLassic implant connections, in different types of transmucosal height.

www.dentaltechitalia.com // info@dental-tech.it
 Facebook: /DentalTechSrl // LinkedIn: /dental-tech-s.r.l.

Surgical Tray

Tray M - For ImpLogic® AT, ImpLassic FT3 and ImpLassic FTP



Dental Tech surgical trays, with their ergonomic design, provide the best ratio between the surgical instrument contained and the size of the box. They are light, manageable, stackable, compact and spacious.

They are also easily transportable and guarantee the safety of their contents, avoiding the loss of the same during handling. The plastic materials used are highly impact resistant and certified for more than 1000 autoclave sterilization cycles.

Tray M Surgical Kit, allows the clinician to work with a single kit of surgical instruments, regardless of the implant line used (ImpLogic® AT, ImpLassic FT3 and ImpLassic FTP).

For easy identification of each instrument and subsequent repositioning after the cleansing and cleaning phase, the professional can count on help of a color-coded system that traces the most suitable surgical procedures for the various implant diameters. All the drills included have practical depth stops that can be easily inserted and removed from the drills in the tip > shank direction. To complete the kit, the screwdrivers and the torque wrench are included. The latter is a high precision tool, easy to use, that provides the professional a tool of tightening and loosening screws, prosthetic elements and implants.

www.dentaltechitalia.com // info@dental-tech.it
 Facebook: /DentalTechSrl // LinkedIn: /dental-tech-s.r.l.