

ImplantBOOK

2023 | Global guide for Dealers and Dentists



DENTAL IMPLANTS p. 5

AROUND IMPLANTOLOGY p. 23


IMPLANTBOOK REPORT p. 55



DENTAL TECH IMPLANT LINES: 1 DESIGN... 3 CONNECTIONS

 **DENTALTECH**



 **DENTALTECH**
INSPIRED BY LITTLE THINGS
inspired by little things

**1 Design
3 Connections**

**SLIM
IMPLANTS**

**ONE PIECE
IMPLANTS**

**SHORT
IMPLANT**

With **over 45 years of experience** in dental implantology area, Dental Tech is an established and constantly growing company in Italy and abroad and specializes in the design, production and distribution of dental implants.

Dental Tech provides professionals with different implant solutions:

- Universal Internal Connection
- Slim Implant
- Tissue Level Implant
- Short Implant
- One-Piece Implant (Bendable or Mini)

In order to satisfy any clinical need, Dental Tech continues its product range's expansion dedicated to the foreign market.

Alongside the implant lines already present, the **conical connection** system will soon be added.

An implant that will have the same morphology as the Tissue Level and Universal Internal Connection implant line, to carry forward the idea that underlies Dental Tech's thinking: **1 design – 3 connections**.

The morphology presents a Spira profile with a Hybrid progress: flat and radiating towards the root, triangular outward; this allows for greater bone penetration in underprepared sites, resulting in bone preservation.

About the connection, the new implant shows one of the most widespread connection geometries in the world:

- Conical sealing & hex connection
- 11° morse Taper locking
- Internal 2.5 hex connection

Thanks to the possibility of using one surgical kit for all 3 connections (compatible with the principal trendy connections in the industry), the work of professionals will be even easier.

For more info:

www.dentaltechitalia.com // info@dental-tech.it

www.dentaltechworldwide.com // info@dentaltechworldwide.com

Facebook: /DentalTechSrl

Linkedin: /dental-tech-s.r.l./dental-tech-worldwide/